

easy-slide

another Arion brand

Easy-Slide® solution

The ultimate donning aid for open toe compression stockings and tights.



Problem definition

For both caregivers and patients, the problem of donning compression stockings and tights, is one, which occurs on a daily basis. For caregivers, donning and removing compression stockings is a major source of physical strain. In addition, patients who are able to carry out the operation independently may find the process an increasingly arduous one. This often has a direct effect on therapy compliance.

Solution

Easy-Slide has been specifically designed to alleviate the problems experienced when donning compression stockings and tights with an open toe. Easy-Slide is made from very smooth material, which has been treated with a special coating. Together with the dual-layer application, this coated material creates a super-smooth aid that reduces friction significantly.

The friction reduction of Easy-Slide helps minimize the physical strain otherwise experienced by caregivers and patients. In turn, this has a positive effect on the therapy compliance of patients. Also, with clear and simple instructions-for-use, patients can enjoy the freedom and convenience of independently donning their patients' compression hosiery for a long time.

Users

Easy-Slide provides an effective solution for patients who wear compression stockings and tights with open toe. Thanks to the Easy-Slide, hosiery can now be donned without difficulty. Furthermore, Easy-Slide is an essential aid for caregivers who are faced with the prospect of donning their patients' compression hosiery on a daily basis.

Easy-Slide®


The attractive alternative

Product benefits


- Reduction of physical strain for both caregiver and user
- Increase in patient's independence and therapy compliance
- Stockings retain therapeutic graduated compression longer
- Easy to use
- Saves time and money
- Durable
- Take it anywhere

Product range

Art.Nr.	size	shoe size
E01500	Extra Small	≤ 34
E01501	Small	34-37
E01502	Medium	38-41
E01503	Large	42-45
E01504	Extra Large	≥ 46



Slide the Easy-Slide over the foot



Slide the stocking over the Easy-Slide



Draw the stocking over the heel



Grasp the loop on the Easy-Slide



Pull the Easy-Slide out